



Innovating Technologies for Life

# DonorCare<sup>®</sup> Needle Guard

Sharps safety for whole blood collection

*Reduces Risk of needlestick injury*

*Immediate Shielding of the needle on withdrawal*

*Facilitates Safe and fast needle withdrawal*

*Improves Safety for subsequent procedures e.g. stripping, segmenting*



[www.itlcorporation.com](http://www.itlcorporation.com)



Innovating Technologies for Life

# DonorCare<sup>®</sup> Needle Guard

## Sharps safety for whole blood collection

DonorCare<sup>®</sup> immediately shields the needle on withdrawal reducing risk of exposure and needlestick injury.


### Safe and Fast

DonorCare<sup>®</sup> allows safe and fast needle withdrawal, even in emergency situations. Fast withdrawal technique increases donor comfort and reduces risk of injury upon withdrawal.

### Easy to Integrate

DonorCare<sup>®</sup> features and procedures are easy to integrate into existing collection processes. DonorCare<sup>®</sup> may be pre-attached on commercially available blood bag sets or simply attached by the end user.

DonorCare<sup>®</sup> Needle Guard is compatible with companion product Samplok<sup>®</sup> Tube Holder. After the whole blood collection is complete, the DonorCare can be inserted into the Samplok<sup>®</sup> where it locks securely providing an added level of safety.



▲ DonorCare<sup>®</sup> is available in two versions to accommodate the various sizes of needle hubs.



▲ DonorCare<sup>®</sup> locked into Samplok<sup>®</sup> Tube Holder provides a safe, secure system prior to disposal in the sharps container.

### Product Information

Latex Free Single Use Only

### Product Codes

**A100100 Small Channel** recommended for use with needle hubs up to 8.5mm wide, 1000 / case

**A100150 Large Channel** recommended for use with needle hubs from 8.5mm - 11mm wide, 1000 / case

DonorCare<sup>®</sup> and SampLok<sup>®</sup> are registered trademarks of ITL Corporation, Melbourne, Australia. Copyright © 2009 ITL Corporation. All rights reserved.

### ITL Corporation

t: +(1) 888 411 2851

f: +(1) 703 435 6717

[www.itlcorporation.com](http://www.itlcorporation.com)  
[sales@itlus.com](mailto:sales@itlus.com)

# Sharps Safety & Best Practice

ITL-BR-0022 REV-1